



## Certificate of Analysis

---

M101: No-Fade™ Mounting Medium

---

**Product:** M101 – No-Fade™ Mounting Medium

**Description:** This fade-retardant medium contains PBS, glycerol, formalin, and 2% DABCO. It is packaged in a 3.5 mL dropper bottle that releases approximately 45 uL per drop, enough medium for 75 applications when applying one drop per well.

**Preparation:** DABCO dissolved in PBS with glycerol, and formalin.

**Lot/Quality Control Number:** 139

**Control/Production Date:** 01/2020

**Stability & Storage:** Store at 4° C.